Chemical Reaction Sensor Lab Name	Per
/as you read; X as you do! This is a lab best done by 2 people. Read out loud. Both wo	ork and mark <u>X</u> for done.
1. Your <u>Problem</u> : Determine if energy is gained or lost during a chemical reaction of NaHCO <sub>3</sub> · You will measure the temperature change during the reaction to find out if energy is gained	+ C <sub>2</sub> H <sub>4</sub> O <sub>2</sub> . or lost.
<ul> <li>A chemical change takes place during a chemical reaction.</li> <li>A chemical change creates a new substance and either uses or releases energy.</li> <li>When sodium bicarbonate (NaCOH<sub>3</sub>) or baking is added to 5% Acetic acid (C<sub>2</sub>H the result is carbon dioxide (CO<sub>2</sub>) gas AND, the solution changes temperature.</li> <li>The chemical equation is NaHCO<sub>3</sub> + C<sub>2</sub>H<sub>4</sub>O<sub>2</sub> → NaCO<sub>3</sub> + H<sub>2</sub>O + CO<sub>2</sub></li> <li>Exothermic reactions give off heat or become w</li> <li>Endothermic reactions take in and use heat so become c</li> </ul>	ł₄O₂ or v)
3. Finish this <u>Hypothesis</u> (as your prediction). Always use a complete sentence with no "it". When vinegar and baking soda react to form carbon dioxide, the t will go showing that energy is The reaction will be _	emperature
(up or down) gained (in) or lost (out)	
4. Be Ready: a) Wear goggles! Follow the directions! b) Protect your ComputerA laptop should be up on something. Cords should be out of sink and any liquid	
Cords should be out of sink and any liquid.     c) Place NO PRESSURE on sensor. Handle sensor with care.     d) Prepare and Test the Application.     Open "Vernier Graphical Analysis".     Chaose Sensor Data Collection.	
<ul> <li>Choose Sensor Data Conection. Prig in your temperature sensor.</li> <li>Use "Mode" to set the Start and Stop on "Manual"</li> <li>Choose "Graph and Table" in the upper right corner box.</li> <li>Each of you use the "Graph Tools" (lower left) to draw your prediction for you</li> <li>"Collect" Take turns placing the sensor inside your elbowWatch "Stop".</li> </ul>	Ir elbow temperature.
<ul> <li>5. Procedures and Observations: Check / as you read. X as you do.</li> <li>a) Open "Graphical Analysis" Start NEW Check ALL of the above settings.</li> <li>b) Sensor should be plugged in.</li> <li>c) Be Safel Goggles should be on No big sleeves Long hair tied back</li> </ul>	
Checked	
<ul> <li>e) Each of you draw your prediction graph on screen. Change "Prediction 1(2)" to your</li> <li>f) Obtain Sodium Bicarbonate on a paper towel. Fold towel to make it easy to pour. DO</li> <li>g) "Collect" to record temperature.</li> <li>b) Record the beginning (before) temperature bare</li> </ul>	Name(s). Save. NOT ADD.
<ul> <li>i) NOW, with one hand holding the flask, <u>SLOWLY</u> add your Sodium bicarbonate (baking soda) to the acetic acid (vinegar).</li> <li>i) Collect until the hubbles are very slow or stop. Click on "Stop"</li> </ul>	g
<ul> <li></li></ul>	Ľ
6. Use Graph Tools ("Edit Graph Options" and "Add Annotations") to label: its name, both prediction names. x-axis & y-axis.	your names,
<ul> <li>Save.</li> <li><u>Conclusion</u>: Add your Conclusion as an Annotation next to the title, using a complete sem When vinegar and baking soda react to form carbon dioxide, th temperature goes showing that energy (in the form of hu</li> </ul>	Checked tence with no "it". re eat)
is showing that it is an reaction.	Save Checked