pH Home Lab Tests	Name	Period
Reference: Central Science→Physical Science→pH Ch	art	
Background Information:		
1) Litmus Paper is used to test for pH:	RED Litmus turns BLUE in B	
	BLUE Litmus turns RED in D.	CHEMISTRY

NO change for both in N

- 2) **pH testing paper** comes in different ranges of pH for different applications of pH testing.

 Our paper tests for pH's from 1 to 12. You must use the chart that matches your testing paper. Use the online chart (Central Science → Physical Science → pH chart)
- 3) pH testing paper requires liquid. Add water for SAFE solid substances.
- 4) Get Permission. BE SAFE. DO NOT USE DANGEROUS SUBSTANCES.



Procedures:

- 1. Obtain pH testing paper and both red and blue Litmus Paper to take home.
- 2. At home, choose safe liquids or solids that you can dissolve in water such baking soda, sugar or soap powder. Fill in the chart as you go.
- 3. PREDICT first!
- 4. Test with a small piece (1-2 cm) of each testing paper, ONE AT A TIME.
- 5. In chart, record FIRST wet COLOR for EACH paper. Record as you go.
- 6 Complete the chart.
- 7. Attach all papers. 2 points each

Substance	Predict	Col	or of	Litmus	Color of	pН	рН	Attach test papers
	:		mus	Results:	pН	number	Paper	here:
		Pa	per		Indicator	from pH	Results:	
	Acid,		n first	Acid,	Paper	Chart	Acid,	
	Base or	W	ıet .	Base or	when		Base or	
	Neutra	0.1		Neutral	first wet		Neutral	
	l5	Rine	Red	,			3	
1								
2								
3								
4								
5								
6								
7								
0								

1	Lave	Cl	-14-4	
	77VP	(.ne	rkon	

Other pH Indicator Liquids: 5+ point	Other	Ha	Indicator	Liauids:	5+	points
--------------------------------------	-------	----	-----------	----------	----	--------

- 1) Choose another juice or safe liquid to use as a pH indicator.
- 2) Make a chart.
- 3) Test some known acids and bases from above, such as vinegar, baking soda & water.
- 4) Record your data.
- 5) If it works, more points for bringing some of your Indicator liquid to class to test.